

## ACE AIROTORQ - Award Winning Technology



**ACE AIROTORQ** is a brushless electric micro-motor, air to electric conversion kit that installs on any delivery system. This most innovative award winning technology from Technodent, redefine modern dental standards with high-precision performance and high-level reliability. The ACE AiroTorq is powered by a Swiss brushless, dual lighting micro-motor technology for seamless results in both high and low speed procedures.

Technodent Exim India Pvt. Ltd.

WEB: [www.technodent.in](http://www.technodent.in) | EMAIL: [cs@technodent.in](mailto:cs@technodent.in) | PHONE: +91 99102 11888 / +91 11 4155 0193

ADDRESS: B138 –GF, NAVAKAR BUILDING MOHAMMADPUR, BHIKAJI CAMA PLACE, NEW DELHI – 110066

The art and science of restorative dentistry is dependent on essential equipment functioning as near to perfection as possible. Practitioners who have tried electric hand-pieces become quick converts to the benefits of preparing a tooth for a crown with an electric hand piece. The torque and versatility of an electric hand-piece are unmatched and allows the practitioner to do better dentistry on his patients.

### Features:

1. Powered by a D-302 – A2E with Blue (UV) light PCB board
2. **Light weight Swiss brushless micro-motor MO-33**
  - a. Offers smooth constant torque (Max 3.5 Ncm) for restorative treatment
  - b. Less noise and vibration
  - c. **Top Feature:** Dual lighting White (25,000 lux) & UV Led for composite cutting
  - d. **Top Feature:** Reduced maintenance requirements with lubrication-free motor
3. External power supply offers high wattage DC 36V 120W for high class performance & safety
4. Included High Speed contra-angle 1:5 hand piece with ceramic bearings & LED in super light titanium sleeve
5. **Top Feature:** Water can be controlled from the hose connector on the motor end
6. Speed control knob / forward & reverse operation
7. **Top Feature:** Installs on the existing delivery system & Operates with the existing Rheostat in the dental unit
8. **Top Feature:** Easy quick connect system for motor cord
9. **Top Feature:** Anti-bacterial & bio-film protected water hose certified by DVGW, Germany
10. Motor speed range from:
  - a. 1500 – 40,000 rpm (1:1)
  - b. 15,000 – 200,000 (1:5)
  - c. 100 – 2500 rpm (16:1)
11. Light, table-top and small control box for easy handling
12. Two years warrantee on motor, tubing & control box
13. After sale service and spares available 24X7

### Installation:





T2 Line is the robust straight and contra-angle hand piece series in the comfort class. Quality and Sirona's expertise are combined to create a product series with functionality and reliability to ensure complete focus on the treatment.

## Optional Instruments:

The Mont Blanc® range of contra-angles delivers all the features you expected!

Recognized for reliability, the Anthogyr contra-angles are user-friendly: the I-Coat® surface treatment delivers optimal tactile sensation in humid environments, its ergonomics delivers perfect handgrip, and the ceramic ball-bearings make use quiet. The combination of two noble materials, titanium and stainless steel, makes it both lightweight and powerful.



## Instruments which can also be used with AiroTorq:

