



EXO SAFE

Automatic periosteal



Composed by a pre-calibrated impaction instrument and blades, the automatic periosteal Exo Safe is a solution for atraumatic dental extractions. Multi-purpose, it becomes an automatic osteotome when used together with the OsteoSafe insert.*



Progressive and controlled impaction

- Maximum bone plate preservation
- Easier mobilization of the tooth and fractured roots
- Effective and time-saving procedure



6 blade types

- Straight convex periosteal wide and thin
- Straight concave periosteal medium and thin
- Angled convex periosteal vertical and horizontal



Universal connection ISO 3944

- Compatible with all types of implantology and chair motors
- One single motor speed 15 000 rpm max ; 1:1 Ratio



Specific tri-lobe connection

- Excellent insert retention during use
- Intuitive and secure locking
- Maximum torque resistance

Products references

COMPLETE EXTRACTION KIT (including 1 impactor Safe + 6 periosteal + 1 sterilization box)	SET OF 6 PERIOTOMES + 1 STERILIZATION BOX
EXO_FULL	EXO_SET

See page 26 for **technical data**.
See page 29 for **accessories**.

*See our "equipments for implantology and surgery" catalogue.