



# IMPULSION® RANGE

Implantology contra-angles



Contra-angles designed for the implant placement: drilling, tapping of the sockets that will receive the dental implants, placement of the dental implants and screwing of the prosthetic parts.



Ref. 14400BP



**Optimised head size\***

→ Good visibility of the implant site for better comfort.



**Easy disassembly using a spanner**

→ Compliance with asepsis standards.  
→ Easier maintenance offering increased longevity.



**Ergonomic design**

→ Excellent handling.



**Depth stop system for controlling the drilling\***

→ Security and comfort while preparing the implant site.



Ref. 14400

Ref. 14401

Product references

	CONTRA-ANGLE LATCH TYPE HEAD	CONTRA-ANGLE PUSH BUTTON HEAD
WITHOUT DEPTH STOP	14400	14400BP
WITH DEPTH STOP	14401	-