



# MONT BLANC® RANGE

General practice equipment



Contra-angles dedicated to dental care (enamel work, elimination of decayed tissues, cavity crowns preparation, removal of filling materials, finishing and polishing of dental and restoration surfaces).



Ref. 10000XL



**Miniaturized head**  
→ Easier access



**Internal quadruple spray irrigation**  
→ Cooling of the site and conservation from heat



**Ceramic bearings\* reducing the wear and tear associated with friction**  
→ Increased strength and durability for enamel work



**Functional design and protective I-coat® coating**  
→ Anatomical handling and scratch resistant



Ref. 10200XL



Ref. 10000XL



Ref. 10300XL

Products references

	CONTRA-ANGLE MULTIPLIER 1:5 (DRILL CA Ø1.6FG)	CONTRA-ANGLE 1:1 (DRILL CA Ø 2.35)	CONTRA-ANGLE REDUCER 5:1 (DRILL Ø 2.35)
LIGHT	10000XL	10200XL	10300XL
WITHOUT-LIGHT	10000X	10200X	10300X

\*Only for ref. 10000X and 10000XL contra-angles.

See page 26 for **technical data**.  
See page 29 for **accessories**.