

RELIABLE - PRACTICAL - ESSENTIAL !

NiTi CONTROL[®] - MICRO NiTi[®]

RANGE OF ENDODONTIC HANDPIECES



*Simple and
precise !*

anthogyr


A global solution for dental implantology

anthogyr's range of endo contra-angles provides practitioners comfort and safety for successful root canal treatment.

The micro head offers maximum visibility and excellent access in the molar area.



Thanks to its adjustment ring, the **NiTi CONTROL®** allows **precise torque control** and thus, avoids files breaking.

Applicable to all standard nickel-titanium file systems.




MAIN ASSETS

- **Durability** thanks to the chrome coating, which improves the resistance of the material.
- **Micro head** for better visibility and better access to the area to be treated.
- **Impeccable maintenance** thanks to an even surface which avoids the retention of bacteria.
- **Quick setting-up of the tool!** Insertion of the file without pressing the push-button.
- **Optimal safety - NiTi CONTROL®**: the automatic declutching function avoids files breaking. Easy and precise adjustment of the torque. Patented system.



Easy adjustment of the torque



	MiCRO NiTi® Micro head push-button	MiCRO NiTi® Latch type head	NiTi CONTROL®	Motor speed required for the files : 300 - 350 rpm	Motor speed required for the files : 600 - 650 rpm
128:1	10628MBP	10628	10628MBP-DA	40 000	
64:1	10664MBP	10664	10664MBP-DA	20 000	40 000
32:1	10632MBP	10632		10 000	20 000
16:1	10616MBP	10616		5 000	10 000
8:1	10608MBP	10608			5 000